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1. 蔡文城, 蔡岳廷。實用臨床微生物診斷學, 第十版。2011: 65-6, 97, 601-5, 653-5。九州圖書文物有限公司, 台北。
2. Roelofsen E, van Leeuwen M, Meijer-Severs GJ *et al.* Evaluation of the effects of storage in two different swab fabrics and under three different transport conditions on recovery of aerobic and anaerobic bacteria. *J Clin Microbiol* 1999; 37:3041-3.
3. Lauer BA, Masters HB. Toxic effect of calcium alginate swabs on *Neisseria gonorrhoeae*. *J Clin Microbiol* 1988; 26:54-6.
4. Mandler F, Sfondrini D. Evaluation of survival of bacteria on dry swabs and transport systems. *Ann Sclavo* 1977; 19:537-45.
5. Rishmawi N, Ghneim R, Kattan R *et al.* Survival of fastidious and nonfastidious aerobic bacteria in three bacterial transport swab systems. *J Clin Microbiol* 2007; 45:1278-83.
6. Perry JL. Assessment of swab transport systems for aerobic and anaerobic organism recovery. *J Clin Microbiol* 1997; 35:1269-71.
7. Clinical and Laboratory Standard Institute (CLSI). Quality control of microbiological transport systems. Approved standard M40-A. 2003. NCCLS, Wayne, PA, USA.
8. Chapel T, Jeffries C, Gatewood C. Effect of medium temperature on recovery of gonococci. *Am J Clin Pathol* 1979; 72:84-7.
9. Ratner HB, Tinsley H, Keller RE *et al.* Comparison of the effect of refrigerated versus room temperature media on the isolation of *Neisseria gonorrhoeae* from genital specimens. *J Clin Microbiol* 1985; 21:127-8.

10. Sng EH, Rajan VS, Yeo KL, Goh AJ. The recovery of *Neisseria gonorrhoeae* from clinical specimens: effects of different temperatures, transport times, and media. *Sex Transm Dis* 1982; 9:74-8.
11. Morosini MI, Loza E, Gutiérrez O *et al.* Evaluation of 4 swab transport systems for the recovery of ATCC and clinical strains with characterized resistance mechanisms. *Diagn Microbiol Infect Dis* 2006; 56:19-24.

Evaluation the Efficacy of CMP™ Aerobic Transtube for Use in the Preservation of Clinical Bacteria

Yen-Hsin Chiu¹, Shian-Shiou Shie¹, Cheng-Chan Lo¹, Yu-June Huang¹,
Shih-Chien Yang¹, Gen-Der Ho¹, Yueh-Ting Tsai^{1,2}

SuperLaboratory Ltd, New Taipei City, Taiwan¹;
Institute of Biochemical Science and Technology, College of Life Science, NTU, Taipei, Taiwan²

Abstract

The efficiency of bacteria preservation of CMP™ Transtube (Creative Biotechnology Co., Ltd. Taiwan.) was evaluated according to the CLSI approved standard M40-A “Quality control of microbiological transport systems”(2003). The experiments used fastidious clinical bacteria, *Neisseria gonorrhoeae* and *Streptococcus pyogenes*, to evaluate the efficacy of CMP™ transtube. After inoculated 10⁶ CFU organisms on the transtube, the transtubes were preserved in room temperature or at 4°C to simulate how the clinical specimen is transported. The survived *N. gonorrhoeae* were recovered by saline after 0, 8, 24 hours storage, and *S. pyogenes* were recovered after 0, 8, 24, 48 hours storage. The recovered bacterial suspensions were performed in a ten-fold serial dilution method, and then plating with spreading method and counting the number of the survival bacteria after incubating at suitable conditions. According to CLSI M40-A statement: “In room temperature, there should be no more than a 3 log decline in CFU between

the zero-time CFU count and the CFU of the swabs that were stored at different time interval; at 4°C, there should be no more than a 1 log increase in CFU and no more than a 3 log decline in CFU between the zero-time CFU and the CFU of the swabs that were stored at different time interval.” The results of this study showed that CMP™ Transtube can well preserved the *S. pyogenes* up to 48 hours whether stored in room temperature or at 4°C. It fulfills the CLSI M40-A guideline. The CMP™ Transtube can also well-preserved *N. gonorrhoeae* up to 24 hours in room temperature, but only last 8 hours at 4°C. The experiment result could be act as reference data when Clinical Microbiology Laboratory choosing their transport device.

Key word: Transport medium, Transtube, *Neisseria gonorrhoeae*, *Streptococcus pyogenes*.